

# BiotinPure™

## Antibody Biotinylation Kit



Easy and fast protocol with no residual biotin

The BiotinPure™ Antibody Biotinylation Kit is a complete kit for easy biotinylation of carrier-free antibodies at room temperature. The biotinylated antibodies are ready for use within one hour.




[www.kromnigon.com](http://www.kromnigon.com)

- Fast protocol  
Biotinylated antibodies in 1 hour
- No residual biotin  
Wash column included
- Everything in one tube  
Reaction inside spin column
- Robust  
Works with most carrier-free antibodies.

# BiotinPure™ Antibody Biotinylation kit

Biotinylate and purify your own carrier-free antibodies in one hour using a standard benchtop centrifuge.




## Antibody Biotinylation Made Easy

Each kit contains enough reagents and consumables to biotinylate five carrier-free antibodies (no protein, glycerol, azide or other chemicals interfering with NHS-ester chemistry in the storage buffer).

BiotinPure™ uses amine/ester chemistry to attach approximately 2 biotins per antibody and includes a purification step that removes unreacted free biotin from the final biotinylated antibody solution.


Biotinylation and purification is completed within one hour using a standard benchtop centrifuge.

## BiotinPure™ Workflow



### BiotinPure™ kits

BiotinPure™ Antibody Biotinylation kit (5 reactions)

 *Developed for use with our StreptaClick® kits, but compatible with many other applications.*

Contact us for custom molecules, bulk orders and OEM opportunities