

Spectrasplit Product List

Product name	Applicable fluorochromes
S-Split DAPI	DAPI. In tissue, this set gives less autofluorescence in ECM.
S-Split 405	DAPI, Dylight350, CF350, AMKA, Brilliant violet 421, AF405, etc.
S-Split 440	Brilliant violet 480, AF440, ATTO440, ATTO425, Sytox blue, etc.
S-Split 488	AF488, ATTO488, CF488, DyLight488, FITC, eGFP, CFSE, etc.
S-Split 546	AF546, ATTO542, Cy3, CF543, AttoRho6G, DyLight549, PE, EFlour570, etc.
S-Split 546/PE	This set activates PE fluorochromes at low energy which prolongs life length.
S-Split 594	AF594, ATTO590, Atto594, Cy3.5, CF594, DyLight594, etc.
S-Split 647	AF647, ATTO647, Cy5.5, DyLight650, etc.
S-Split 647/APC	Brighter and shifted some to work excellent for both AF647 & APC.
S-Split 750	Fluorochromes here usually requires a LED-lamp.

