

BiotinPure™

Easy antibody biotinylation

Fast protocol - no residual biotin

The BiotinPure™ Antibody Biotinylation Kit is a complete kit for easy biotinylation of antibodies at room temperature within one hour.



www.kromnigon.com



- Fast protocol**
Biotinylated antibodies in 1 hour
- No residual biotin**
Wash column included
- Everything in one tube**
Reaction inside spin column
- Robust**
Works with any buffer without stabilizing proteins

BiotinPure™ Antibody Biotinylation kit

Biotinylate and purify your own antibodies in one hour




Antibody Biotinylation Made Easy

Each kit contains enough reagents and consumables to biotinylate five antibodies. BiotinPure™ is optimized for our StreptaClick® reagents, but the biotinylated antibodies can be used for other applications.

BiotinPure™ uses amine/ester chemistry to attach approximately 2 biotins per antibody and includes a purification step that removes unreacted free biotin from the final biotinylated antibody solution.

Biotinylation and purification is completed within one hour.

BiotinPure™ Workflow



BiotinPure™ kits

BiotinPure™ Antibody Biotinylation kit (5 reactions)

Contact us for custom molecules,
bulk orders and OEM opportunities